



SCIENTIFIC COMPUTING AND IMAGING INSTITUTE



ShapeWorks Documentation

ShapeToolKit : Alignment Tools Documentation

Shireen Elhabian, Praful Agrawal, Riddhish Bhalodia,
Joshua Cates, Manasi Datar, Ross Whitaker

July 22, 2018

Abstract

This document covers the alignment and registration based tools for both images and meshes.

Contents

1	License	3
2	Acknowledgement	3
3	Alignment Tools	3
3.1	Center of mass alignment (images)	3
3.2	Rigid ICP based 3D image registration	4
3.3	Rigid ICP based surface mesh registration	4

1 License

ShapeWorks and its graphical user interface ShapeWorksStudio are available for free and are open source under the MIT License. For more information, please refer to <https://www.sci.utah.edu/software/shapeworks.html>.

The MIT License Copyright (c) 2012 Scientific Computing and Imaging Institute, University of Utah. License for the specific language governing rights and limitations under Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions: The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software. THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

2 Acknowledgement

If you use ShapeWorks and/or ShapeWorksStudio in work that leads to published research, we humbly ask that you add the following to the 'Acknowledgments' section of your paper: *"This project was supported by the National Institute of General Medical Sciences of the National Institutes of Health under grant number P41 GM103545-18."*

3 Alignment Tools

3.1 Center of mass alignment (images)

Takes in a binary volume and translates the center of mass of the shape to the center of the 3D volume space or to any user defined center.

The command line tool can be used as follows:

```
./TranslateShapeToImageOrigin -inFilename $1 -outFilename $2  
-MRIinFilename $3 -MRIoutFilename $4 -useCenterOfMass $6  
-centerX $7 -centerY $8 -centerZ $9
```

The parameters are as follows

1. inFilename : Input image file name
2. outFilename : Output image file name
3. MRIinFilename : MRI image file name
4. MRIoutFilename : MRI output image file name
5. useCenterOfMass : Flag to be set as 1 to use center of mass, if 0 then it translates the image center
6. centerX : (Optional) Specified X coordinate of center of the final output (Default = Image Center)
7. centerY : (Optional) Specified Y coordinate of center of the final output (Default = Image Center)
8. centerZ : (Optional) Specified Z coordinate of center of the final output (Default = Image Center)

3.2 Rigid ICP based 3D image registration

```
./ICPRigid3DImageRegistration -targetDistanceMap $1  
-sourceDistanceMap $2 -sourceSegmentation $3 -icpIterations $4  
-isoValue $5
```

A command line tool that performs iterative closed point (ICP) 3D rigid registration on a pair of images.

It uses the following input arguments:

1. targetDistanceMap: the distance map of target image.
2. sourceDistanceMap: the distance map of source image.
3. sourceSegmentation: the segmentation of source image.
4. isoValue: as we need to get point set from surface for ICP, this iso value is required to get the isosurface. The default value is 0.0.
5. icpIterations: the number of iterations user want to run.

We in the provided automated scripts use ANTS based packages to perform registration.

3.3 Rigid ICP based surface mesh registration

Performs iterative closed point (ICP) rigid registration on a pair of vtk meshes.

It uses a parameter file that would enable to specify the source mesh (moving) and the target mesh (fixed) to be used to estimated the rigid transformation matrix then apply the same transformation on other meshes defined in the source mesh domain to be mapped to the target domain. The usage therefore is just

`./ICPRigid3DMeshRegistration parameterFile.xml` parameter file tags are as follows:

1. `source_mesh`: vtk filename of the moving mesh
2. `target_mesh`: vtk filename of the fixed mesh
3. `out_mesh` : vtk filename of the aligned moving mesh to be save
4. `out_transform` : txt filename to save the estimated transformation
5. `source_meshes`: (optional) a list of vtk filenames for meshes defined in the source mesh domain to be mapped to the target domain using the same transformation matrix estimated.
6. `out_meshes` : a list vtk filenames to save source_meshes after applying the transformation matrix.
7. `mode` : Registration mode rigid, similarity, affine (default: similarity)
8. `icp_iterations`: number of iterations
9. `debug`: verbose debugging information
10. `visualize`: display the resulting alignment